


# Size Selection Guide - Which bath should I choose?



## Key

### Unstirred Bath

**8 14 28 56**

### Compact Unstirred Bath

**28 56**

### Duobath™

**8 8 14 14**

### Compact Duobath™

**8 28 8 28**

### Immersion Circulator & Carry Case



## Tips -

**Flexibility** ~ Use Duobath™ to cook at two temperatures simultaneously e.g. fish at 48°C and regen at 60°C.

**Capacity** ~ This may mean the ability to cook larger mis en place e.g. pork belly or the ability to cater for larger number of covers e.g. special event

**Limited Space** ~ Use immersion circulator then remove and store away when not in use, or consider the compact versions FLE28D & FLE56D (see website and brochure for details.)

Stirred baths also available in 28 and 56 litres. Complete with rotating controller and flat hinged bi-fold lid.